

Speckle Instrument GUI – Programmers Guide

Dave Mills (rfactory@theriver.com) – July 2018

Table of Contents

1. Development Environment.....	1
1.1 Dependencies.....	1
1.2 Hardware.....	2
1.3 Serial Numbers.....	3
2. Software Architecture.....	4
3. Directory structure.....	7
4. Tcl wrappers.....	8
5 Code documentation.....	9

1. Development Environment

The recommended development environment is Linux x86_64 with a 3.x series kernel and the normal development toolchain installed. Initial development took place on a CentOS 7 system , but work was also done on Ubuntu 18. Any modern Linux variant should suffice.

1.1 Dependencies

The following package should be installed and configured before working on the code.

- automake
- autogen
- cfitsio*
- doxygen
- ds9
- fftw3*
- gcc
- g++
- glib-2.0*
- gmodule-2.0*
- gnuplot
- gobject-2.0*
- libpng*
- make*
- mysql-client
- mysql-server
- mysqltcl
- pangoft2*
- pkg-config

```
python
tcl*
tk*
xpa-tools
zlib*
```

and their associated -dev or -devel packages for the ones marked with *.

The binaries and libraries for ds9, xpa-tools, tcl, tk are included in the Speckle code distribution, so those can be used as-is with a fully compatible system.


The VIPS source tree is also included and will almost certainly need recompiling in most cases.

1.2 Hardware

A minimum of 2 USB 3.0 capable ports are required to interface the cameras.

Another 3 USB 2.0 or better ports are needed to interface the Filter Wheels and the Zaber motion control stages. All the Zaber devices are daisy-chained off a single port.

Once the system is configured , use lsusb to examine the device complement



```
Terminal
File Edit View Terminal Tabs Help
nessi:lsusb
Bus 001 Device 002: ID 8087:8001 Intel Corp.
Bus 001 Device 001: ID 1d6b:0002 Linux Foundation 2.0 root hub
Bus 003 Device 016: ID 05e3:0612 Genesys Logic, Inc. Hub
Bus 003 Device 017: ID 136e:0012 Andor Technology Ltd.
Bus 003 Device 018: ID 136e:0012 Andor Technology Ltd.
Bus 003 Device 001: ID 1d6b:0003 Linux Foundation 3.0 root hub
Bus 002 Device 003: ID 8087:0a2a Intel Corp.
Bus 002 Device 044: ID 104d:1011 Newport Corporation
Bus 002 Device 043: ID 104d:1011 Newport Corporation
Bus 002 Device 042: ID 0403:6001 Future Technology Devices International, Ltd FT
232 Serial (UART) IC
Bus 002 Device 041: ID 05e3:0610 Genesys Logic, Inc. 4-port hub
Bus 002 Device 001: ID 1d6b:0002 Linux Foundation 2.0 root hub
nessi:
nessi:
nessi:
nessi:
nessi:
nessi:
nessi:
nessi:
nessi:
nessi:
nessi:
```

The Filter Wheels show up as Newport Corporation, and the Zabers as Future Technology....
In order to control these devices as a non-root user run the script

./setDevicePermissions

in the base folder of the Speckle software installation.

1.3 Serial Numbers

If it becomes necessary to change out either Filter Wheel or Camera components, the appropriate configuration files will be adjustment. The configuration files are in the \$HOME/speckle-control directory

andorsConfiguration.[telescope]

filtersConfiguration.[telescope]

In each case the serial number information will need to be updated.
The Filter Wheel serial numbers can be found using the lsusb command

```
nessi:lsusb
Bus 001 Device 002: ID 8087:8001 Intel Corp.
Bus 001 Device 001: ID 1d6b:0002 Linux Foundation 2.0 root hub
Bus 003 Device 016: ID 05e3:0612 Genesys Logic, Inc. Hub
Bus 003 Device 017: ID 136e:0012 Andor Technology Ltd.
Bus 003 Device 018: ID 136e:0012 Andor Technology Ltd.
Bus 003 Device 001: ID 1d6b:0003 Linux Foundation 3.0 root hub
Bus 002 Device 003: ID 8087:0a2a Intel Corp.
Bus 002 Device 044: ID 104d:1011 Newport Corporation
Bus 002 Device 043: ID 104d:1011 Newport Corporation
Bus 002 Device 042: ID 0403:6001 Future Technology Devices International, Ltd FT
232 Serial (UART) IC
Bus 002 Device 041: ID 05e3:0610 Genesys Logic, Inc. 4-port hub
Bus 002 Device 001: ID 1d6b:0002 Linux Foundation 2.0 root hub
nessi:
nessi:
nessi:lsusb -v -s 002:043 | grep iSerial
    iSerial          128 061D088E010F5400
nessi:lsusb -v -s 002:044 | grep iSerial
    iSerial          128 1B18177A01135400
nessi:
```

The Andor Serial numbers can be found by examining the “dmesg” log at system boot time.

2. Software Architecture

The Speckle GUI and control package includes the following components

- Tcl/Tk scripts for GUI and controls
- Shared Libraries for Andor and Filter Wheel controlled
- Shared libraries for FITS I/O and Image processing
- ds9 and XPA for Image display

When the software is running there are normally 6 processes involved

- ds9red - displays the red camera images
- ds9blue - displays the blue camera images
- andorCameraServer 0 - controls Andor camera with id=0
- andorCameraServer 1 - controls Andor camera with id =1
- gui2 - GUI and control windows
- xpans - XPA protocol server

The following communication methods are used

- The GUI interacts with the camera servers via sockets (2001 and 2002)
- It is also possible to telnet to these sockets and send camera server commands manually.

- The GUI interacts with the ds9 displays via XPA (using xpaget, xpaset programs)

- The Camera servers interact with the ds9 displays via shared memory buffers

- The GUI interacts with the Camera servers via a shared memory area
(a small number of items are used to communicate during the high speed acquisition loop)

Some prototype code is included in the Filter Wheel and Zaber scripts to facilitate their use in socket based server mode (Not implemented).

Separating the Camera low level control into servers make it possible to reset a misbehaving camera without affecting the rest of the system.

The shared memory areas can be listed using the ipcs tool

```
Terminal
File Edit View Terminal Tabs Help
nessi:ipcs -a

----- Message Queues -----
key          msqid      owner      perms      used-bytes   messages


----- Shared Memory Segments -----
key          shmid      owner      perms      bytes        nattch     status
0x00000000   65536      nessi      600        393216        2          dest
0x00000000   98305      nessi      600        393216        2          dest
0x00000000   327682     nessi      600        393216        2          dest
0x00000000   79298563   nessi      600        524288        2          dest
0x00000000   622596     nessi      600        524288        2          dest
0x00000000   753669     nessi      600        524288        2          dest
0x00000000   786438     nessi      600        393216        2          dest
0x00000000   917511     nessi      600        524288        2          dest
0x00000000   950280     nessi      600        524288        2          dest
0x00000000   1212425    nessi      600        12288         2          dest
0x00000000   1245194    nessi      600        393216        2          dest
0x00000000   1277963    nessi      600        12288         2          dest
0x00000000   1310732    nessi      600        393216        2          dest
0x00000000   1343501    nessi      600        12288         2          dest
0x00000000   1376270    nessi      600        393216        2          dest
0x00000000   1409039    nessi      600        12288         2          dest
0x0000d9b5   78741520   nessi      666         488           3
0x00001e5c   1474577    nessi      666       4194304        2
0x00000000   1605650    nessi      600        524288        2          dest
0x00001e5b   1638419    nessi      666       4194304        2
0x00000000   6520852    nessi      600        393216        2          dest
0x00000000   2588693    nessi      600        393216        2          dest
0x00000000   10289174   nessi      600        393216        2          dest
0x00000000   10387479   nessi      600        393216        2          dest
0x00000000   10420248   nessi      600        524288        2          dest
0x00001e5d   19726361   nessi      666         56            3

----- Semaphore Arrays -----
key          semid      owner      perms      nsems

nessi:
```

For Speckle, the ones involved have keys 1e5b, 1e5c, and 1e5d, and each will have an nattch count of 2 during normal operations.

The 2 large areas contain the image buffers, each has sufficient space for a full-frame image at 4 bytes per pixel (1024 x 1024 x4).



The 3rd small area contains a number of control and reporting variables.

```

typedef struct shmControlRegisters {
    int iPeak[2];
    int iMin[2];
    int iFrame[2];
    int displayFFT;
    int displayLucky;
    int saveLucky;
    int iabort;
    int iLuckyThresh[2];
    int iLuckyCount[2];
} shmControl;

```

3. Directory structure

The base directory is

speckle-control

It holds the configuration files, startup scripts, and gemini library archive.

The subdirectories are

- andor – the Andor camera scripting and C code
- bin – executables
- ccd – Utility ccd data processing library and wrappers
- doc – Device and code documentation
- gui-scripts – GUI and control scripts
- include – Include files for re-compilation of libraries
- lib – shared libraries and tcl packages
- oriel – Filter wheel control scripts and library
- picomotor – Pico motor control scripts
- share – installed by VIPS packaged
- vips-8.5.9 – VIPS sources

The contents of most of the subdirectories is self explanatory.

The andor subdirectory contains the following

- andorCameraServe.tcl – The main camera server script (2 instances of this are run)
- andorCodeGen.tcl – Script to auto-generate wrappers for most Andor library calls (Set/Get)
 - Produces
 - andorCreateTclCmds.c
 - andorGenTclInterfaces.c
 - andorGenTclInterfaces.h
- andor.tcl – Scripts to facilitate communication with GUI
- andor_tcl.c – Wrappers for data acquisition and processing, andor functions
- andor_tcl.h - Wrappers for data acquisition and processing, andor functions
- atcmdLXd.h – Andor library function definitions
- buildandorWrap – script to rebuild andorTclInit.so
- ccd_astro.c – Utility astronomy functions
- dofft.cpp – Uses VIPS to perform FFT on image data-directory
- ds9refersher.tcl – Loaded into ds9's to perform fast image refresh
- examples – Andor provided example code

4. Tcl wrappers

All the low-level C and C++ software is accessed from the tcl scripting layer by means of wrappers. The wrappers are C code which is compiled with the tcl api and the resulting library is loaded at runtime into the wish executable (Tcl/Tk shell program).

In order to add a new command , the following steps need to be performed

e.g. To add a command to the andor interface (compiled into andorTclInit.so)

1. Edit the file andor/andor_tcl.c to add a prototype

```
/**
 * \brief tcl_andorMyNewCommand A very useful addition
 * \param ClientData Tcl handle
 * \param Tcl_Interp interpreter pointer
 * \param argc Argument count
 * \param argv Arguments
 */
int tcl_andorMyNewCommand(ClientData clientData, Tcl_Interp *interp, int argc, char **argv);
```

2. Inside the Andortclinit_Init function body, add a line to define the new tcl command

```
Tcl_CreateCommand(interp, "andorMyNewCommand", (Tcl_CmdProc *) tcl_andorMyNewCommand, NULL, NULL);
```

3. Write the code for tcl_andorMyNewCommand and add it to the file

```
/**
 * Parameters are passed by Tcl and must be decoded from strings
 * \param width Width of area
 * \param height Height of area
 */
int tcl_andorMyNewCommand(ClientData clientData, Tcl_Interp *interp, int argc, char **argv)
{
    int width,height;

    if (argc < 3) {
        Tcl_AppendResult(interp, "wrong # args: should be \"",argv[0],\" width height\"", (char *)NULL);
        return TCL_ERROR;
    }

    sscanf(argv[1],"%d",&width);
    sscanf(argv[2],"%d",&height);
    printf "Got %d %d",width,height");
    return TCL_OK;
}
```


4. Compile and rebuild the library following the steps in buildandorWrap

```
gcc -g -c -fPIC andor_tcl.c -fpic -DLINUX -DWall -g -I./include -I./include/tcl $(pkg-config --cflags vips)
g++ -g -shared -o andorTclInit.so ccd_astro.o andor_tcl.o andorGenTclInterfaces.o andorCreateTclCmds.o \
    dofft.o -L../lib -lcfitsio ../lib/libandor.so.2 ../lib/libUSB12C.so.2 $(pkg-config --libs vips)
mv andorTclInit.so ../lib/.
```

5. Test

```
wish
load $env(SPECKLE_DIR)/lib/andorTclInit.so
andorMyNewCommand 123 456
```

6. Rebuild the code documentation

```
cd $SPECKLE_DIR/doc/code
doxygen -g specklecode
```

5 Code documentation

Detailed code documentation showing all the routines and relationships between them can be found in

`$SPECKLE_DIR/doc/code/html/index.html`