

Speckle Instrument GUI - Linux User Guide

Dave Mills (rfactory@theriver.com) – June 2018

Contents

1. Introduction
2. Installation
3. Graphical User Interface
 - The Main window
 - Normal usage
 - Engineering mode
4. Remote control
5. External packages

1. Introduction

The Speckle Instrument GUI has been developed by The Random Factory (Tucson, AZ) in collaboration with the Speckle Instrument PI (Steve Howell) and collaborators (Nic Scott, and Mark Everett - KPNO).

2. Installation

The GUI and accompanying packages are packaged using the **gzipped tar** archives. To install the package :

- `tar xvzf speckle-control-x.y.z.tgz`

where x.y.z is the appropriate version number.


This installation will place the files in the directory `./speckle-control`. Although it is possible to install the software to a different location, this is not recommended as it will be necessary to manually change the location in some of the scripts included with the drivers.

Run the Andor drivers installation script

`./speckle-control/andor/TBD`

Configure the USB devices for rw accessed

`./speckle-control/setDevicePermissions`



```
Terminal
File Edit View Terminal Tabs Help
nessi:lsusb
Bus 001 Device 002: ID 8087:8001 Intel Corp.
Bus 001 Device 001: ID 1d6b:0002 Linux Foundation 2.0 root hub
Bus 003 Device 016: ID 05e3:0612 Genesys Logic, Inc. Hub
Bus 003 Device 017: ID 136e:0012 Andor Technology Ltd.
Bus 003 Device 018: ID 136e:0012 Andor Technology Ltd.
Bus 003 Device 001: ID 1d6b:0003 Linux Foundation 3.0 root hub
Bus 002 Device 003: ID 8087:0a2a Intel Corp.
Bus 002 Device 044: ID 104d:1011 Newport Corporation
Bus 002 Device 043: ID 104d:1011 Newport Corporation
Bus 002 Device 042: ID 0403:6001 Future Technology Devices International, Ltd FT
232 Serial (UART) IC
Bus 002 Device 041: ID 05e3:0610 Genesys Logic, Inc. 4-port hub
Bus 002 Device 001: ID 1d6b:0002 Linux Foundation 2.0 root hub
nessi:
nessi:
nessi:
nessi:
nessi:
nessi:
nessi:
nessi:
nessi:
nessi:
```

Once this setup has been completed, the interface can be started with the command

~/startspeckle2

These last two steps can also be done by double-clicking the appropriate desktop icons.

3. Graphical user interface.

The graphical user interface provides easy access to the major functions such as image acquisition, temperature control, and device setup and configuration.


The program may be started from an xterm by typing (in the speckle-control directory)

./startspeckle2

The program will open a small main window, and then create a message window which shows the progress of the system startup operations.

Once the message window closes, the system is ready for use. The camera is initialized, and temperature control has been switched on.

3.1 The main window



AcquireBright
Standard
SetupSingle
SpeckleVideo
SetupSeries
AcquireFaint
User selected
Save current as

Cooler on
Cooler off
Cooler to ambient



Acq-roi-128
Acq-roi-256
Acq-roi-512
Acq-full
Adjust ROI
Reset full-frame

Engineering
Observing
Filter Selection
Camera status
Plot timings
HOME all stages
zabers to wide mode
zabers to speckle mode
zaber red wide
zaber red speckle
zaber blue wide
zaber blue speckle
zaber input wide
zaber input speckle

Camera Configurations

		Blue Arm	Red
Shutter	0	0	
FrameTransferMode	1	1	
OutputAmplifier	0	0	
EMAdvanced	1	1	
EMCCDGain	2	2	
HSSpeed	1	1	
VSSpeed	1	1	
PreAmpGain	1	1	
ReadMode	4	4	
AcquisitionMode	1	1	
KineticCycleTime	0.0	0.0	
NumberAccumulations	1	1	
NumberKinetics	1	1	
AccumulationCycleTime	0.0	0.0	
TExposure	0.060000	0.060000	
TAccumulate	1.267853	1.267853	
TKinetics	1.267853	1.267853	

Close



Speckle Control

File

Set ROI's

Temperature

Tools

Configurations

Help

Exposure0.06Binning1

Num. Frames1Display FFT

Exp. TypeObjectKinetic mode

Num. Seq.1Accum.1

File name : N20180706_0000016

ObserveVideoAbort

Simulate : AndorsZabersTelemetryFilters

Configure data directory

-52.9 degC-54.3 degC

☒ BLUE ARM☒ RED ARM

Comments : ☐ Auto-clear

targetHR 7117

ProgIDtest

ra+09:00:11.469

dec+82:55:17.328

telescopeWIYN

instrumentspeckle

DQ - image

DQ - cloud

DQ - water

DQ - bg

BLUE ARM☒ Frame Transfer

RED ARM☒ Frame Transfer

Shutter=autoTemp Set-60Autofit ds9

Shutter=autoTemp Set-60Autofit ds9

Mode=wideFilter = clear

Mode=wideFilter = clear

☐ EMCCD☐ High GainEM Gain0Auto Set

☐ EMCCD☐ High GainEM Gain0Auto Set

Observation status : Idle

Vspeed1.13

EMCCD HS20

CCD HS1 MHz

Vspeed1.13

EMCCD HS20

CCD HS1 MHz

INPUT


Move to0Move to0Move to0



Set WIDE to currentSet WIDE to currentSet WIDE to current


Set SPECKLE to currentSet SPECKLE to currentSet SPECKLE to current

Set HOME to currentSet HOME to currentSet HOME to current

LoadSaveLoadSaveLoad







A terminal window titled "Terminal" with a menu bar containing "File", "Edit", "View", "Terminal", "Tabs", and "Help". The terminal displays the command "nessi:more /tmp/timings" followed by a list of 31 numerical values, each on a new line. The values range from 0.0 to 2.193000078201294. The terminal has a standard window title bar with minimize, maximize, and close buttons, and a vertical scrollbar on the right side.

```
nessi:more /tmp/timings
0.0
0.06500005722045898
0.12899994850158691
0.1940000057220459
0.25800013542175293
0.3229999542236328
0.3880000114440918
0.4519999027252197
0.5160000324249268
0.5810000896453857
0.6459999084472656
0.7100000381469727
0.7739999294281006
0.8389999866485596
0.9030001163482666
0.9679999351501465
1.0329999923706055
1.0970001220703125
1.1610000133514404
1.2260000705718994
1.2899999618530273
1.3550000190734863
1.4189999103546143
1.4839999675750732
1.5480000972747803
1.6129999160766602
1.6770000457763672
1.7409999370574951
1.805999994277954
1.871000051498413
1.934999942779541
2.0
2.064000129699707
2.128000020980835
2.193000078201294
```

```
nessi:ls /tmp/speckl*.log
/tmp/speckleLog_1530322241.log /tmp/speckleLog_1530502749.log
/tmp/speckleLog_1530322463.log /tmp/speckleLog_1530502891.log
/tmp/speckleLog_1530322723.log /tmp/speckleLog_1530503583.log
/tmp/speckleLog_1530323111.log /tmp/speckleLog_1530504217.log
/tmp/speckleLog_1530323495.log /tmp/speckleLog_1530504290.log
/tmp/speckleLog_1530323737.log /tmp/speckleLog_1530504332.log
/tmp/speckleLog_1530324010.log /tmp/speckleLog_1530504425.log
/tmp/speckleLog_1530324888.log /tmp/speckleLog_1530504611.log
/tmp/speckleLog_1530325368.log /tmp/speckleLog_1530504737.log
/tmp/speckleLog_1530325752.log /tmp/speckleLog_1530504795.log
/tmp/speckleLog_1530326089.log /tmp/speckleLog_1530505286.log
/tmp/speckleLog_1530327038.log /tmp/speckleLog_1530505602.log
/tmp/speckleLog_1530327151.log /tmp/speckleLog_1530506620.log
/tmp/speckleLog_1530332355.log /tmp/speckleLog_1530508085.log
/tmp/speckleLog_1530335084.log /tmp/speckleLog_1530508483.log
/tmp/speckleLog_1530335813.log /tmp/speckleLog_1530511001.log
/tmp/speckleLog_1530336849.log /tmp/speckleLog_1530511098.log
/tmp/speckleLog_1530398480.log /tmp/speckleLog_1530511224.log
/tmp/speckleLog_1530398755.log /tmp/speckleLog_1530546285.log
/tmp/speckleLog_1530398794.log /tmp/speckleLog_1530546434.log
/tmp/speckleLog_1530398814.log /tmp/speckleLog_1530548907.log
/tmp/speckleLog_1530399172.log /tmp/speckleLog_1530553512.log
/tmp/speckleLog_1530399516.log /tmp/speckleLog_1530553964.log
/tmp/speckleLog_1530403411.log /tmp/speckleLog_1530554786.log
/tmp/speckleLog_1530403900.log /tmp/speckleLog_1530556213.log
/tmp/speckleLog_1530404046.log /tmp/speckleLog_1530556852.log
/tmp/speckleLog_1530405721.log /tmp/speckleLog_1530557089.log
/tmp/speckleLog_1530407045.log /tmp/speckleLog_1530557647.log
/tmp/speckleLog_1530410430.log /tmp/speckleLog_1530560460.log
/tmp/speckleLog_1530410786.log /tmp/speckleLog_1530561133.log
/tmp/speckleLog_1530410917.log /tmp/speckleLog_1530561509.log
/tmp/speckleLog_1530411091.log /tmp/speckleLog_1530562077.log
```


4. Command Line usage

There is a rich set of commands to allow interactive and scripted usage.

To access the command line it is necessary to source the *startspeckle-cmds* script from the speckle-control directory and then type

```
source gui-scripts/gui2.tcl
```

The following commands are available

Filter Wheel :

```
loadFiltersConfig [filename]  
saveFiltersConfig [filename]  
echoFiltersConfig  
selectfilter arm filter-number  
findWheels  
resetFilterWheel arm
```

Zaber stages :

```
loadZaberConfig [filename]  
saveZaberConfig [filename]  
echoZaberConfig  
zaberPrintProperties  
zaberConnect  
zaberDisconnect  
homeZabers  
zaberCheck
```

zaberSetPos name position
zabersStopAll
zaberGoto name station

Pico Stages :

loadPicosConfig [filename]
savePicosConfig [filename]
echoPicosConfig
picosConnect
picoCommand axis cmd
picoSet axis parameter value

Andor Cameras :

Command may be issued from the GUI command line, scripted , or optionally by telnet to ports 2001, 2002. When using the command line the syntax is

commandAndor arm “command and parameters”
or commandAndors “command and parameters”

accumulationcycletime seconds
acquisition index
autofitds9 0/1
baseclamp 0/1
comments comment1|comment2|...
configure hbin vbin vstart vend hstart hend preamp vsspeed ccdhss emccdhss
datadir data-directory
dqtelemetry rawiq rawcc raqwv rawbg
emadvanced index
emccdgain 0/1
fastVideo exposure xs ys dim

fitsbits data-format
forceroi xs xe ys ye
frametransfer index
gettemp
grabcube exposure xs ys dim
grabframe exposure
grabroi exposure xs ys dim
hsspeed amp index
imagename image-name
kineticcycletime seconds
locatestar smooth dim
numberaccumulations count
numberkinetics count
outputamp index
positiontelem input-zaber field-zaber filter
preampgain index
programid program-id
readmode index
reset mode
setexposure seconds
setframe mode
setroi mode
settemperature degrees
shutdown
shutter index
status
version
vsamplitude index
vsspeed index
whicharm

5. Recompiling the shared libraries

Low level functionality is provided in C/C++ for speed , and this code is

wrapped using tcl and loaded into the interpreter at runtime.

To move the code to a different version of Linux it may be necessary to recompile the libraries in the following directories. Each has either a Makefile or a set of build steps (e.g. andor/buildAndorWrap).

The Vips library may present more difficulty due to its many dependencies.

The package can be recompiled using the GNU standard incantations

```
./configure --prefix=some-installation-directory --without-python  
make install
```

If the configure step does not work , try

```
sudo apt install automake autogen m4  
libtoolize  
aclocal  
automake --add-missing  
autoconf
```

then try the ./configure step again.

```

nessi:lsusb
Bus 001 Device 002: ID 8087:8001 Intel Corp.
Bus 001 Device 001: ID 1d6b:0002 Linux Foundation 2.0 root hub
Bus 003 Device 016: ID 05e3:0612 Genesys Logic, Inc. Hub
Bus 003 Device 017: ID 136e:0012 Andor Technology Ltd.
Bus 003 Device 018: ID 136e:0012 Andor Technology Ltd.
Bus 003 Device 001: ID 1d6b:0003 Linux Foundation 3.0 root hub
Bus 002 Device 003: ID 8087:0a2a Intel Corp.
Bus 002 Device 044: ID 104d:1011 Newport Corporation
Bus 002 Device 043: ID 104d:1011 Newport Corporation
Bus 002 Device 042: ID 0403:6001 Future Technology Devices International, Ltd FT
232 Serial (UART) IC
Bus 002 Device 041: ID 05e3:0610 Genesys Logic, Inc. 4-port hub
Bus 002 Device 001: ID 1d6b:0002 Linux Foundation 2.0 root hub
nessi:
nessi:
nessi:lsusb -v -s 002:043 | grep iSerial
    iSerial          128 061D088E010F5400
nessi:lsusb -v -s 002:044 | grep iSerial
    iSerial          128 1B18177A01135400
nessi:
0x00000000 1277963   nessi      600      12288      2      dest
0x00000000 1310732   nessi      600      393216     2      dest
0x00000000 1343501   nessi      600      12288      2      dest
0x00000000 1376270   nessi      600      393216     2      dest
0x00000000 1409039   nessi      600      12288      2      dest
0x0000d9b5 78741520  nessi      666       488        3
0x00001e5c 1474577   nessi      666     4194304     2
0x00000000 1605650   nessi      600     524288      2      dest
0x00001e5b 1638419   nessi      666     4194304     2
0x00000000 6520852   nessi      600     393216      2      dest
0x00000000 2588693   nessi      600     393216      2      dest
0x00000000 10289174  nessi      600     393216      2      dest
0x00000000 10387479  nessi      600     393216      2      dest
0x00000000 10420248  nessi      600     524288      2      dest
0x00001e5d 19726361  nessi      666       56         3

----- Semaphore Arrays -----
key          semid      owner      perms      nsems
nessi:

```

TOPCAT

File Views Graphics Joins Windows VO Interop Help

Table List
1: select * from Speckle

Current Table Properties
Label: select * from Speckle_Observations
Location: select * from Speckle_Observations
Name: jdbc:mysql://localhost/speckle#select * from Speckle_Observations
Rows: 5,818
Columns: 25
Sort Order:
Row Subset: All
Activation Action: (no action) ☐ Broadcast Row

SAMP
Messages: Clients:

TOPCAT(1): Table Browser

Window Subsets Help

Table Browser for 1: select * from Speckle_Observations

e	exposureStart	exposureEnd	Filter	amplifier	numex...	numAc...	window	colBin...	rowBin...	RA	Declination	dqImage	dqCloud	dqWater	dqBG
5800	1.530573E9	1.530573E9	clear	ECMMD Amplifier	1	1	1,1024,1,1024	1	1	+20:53:02.780	+87:52:38.838	0	0	0	0
5801	1.530573E9	1.530573E9	clear	ECMMD Amplifier	1	1	1,1024,1,1024	1	1	+20:54:14.16	+87:52:39.23	0	0	0	0
5802	1.530573E9	1.530573E9	clear	ECMMD Amplifier	1	1	1,1024,1,1024	1	1	+20:54:14.16	+87:52:39.23	0	0	0	0
5803	1.530578E9	1.530578E9	clear	ECMMD Amplifier	1	1	1,1024,1,1024	1	1	+21:44:54.319	+88:24:28.219	0	0	0	0
5804	1.530578E9	1.530578E9	clear	ECMMD Amplifier	1	1	1,1024,1,1024	1	1	+21:44:54.319	+88:24:28.219	0	0	0	0
5805	1.530915E9	1.530915E9	clear	ECMMD Amplifier	1	1	1,1024,1,1024	1	1	+10:13:13.532	+82:56:18.104	0	0	0	0
5806	1.530915E9	1.530915E9	clear	ECMMD Amplifier	1	1	1,1024,1,1024	1	1	+10:13:17.85	+82:56:18.141	0	0	0	0
5807	1.530917E9	1.530917E9	clear	ECMMD Amplifier	1	1	1,1024,1,1024	1	1	+10:41:19.991	+82:56:33.363	0	0	0	0
5808	1.530917E9	1.530917E9	clear	ECMMD Amplifier	1	1	1,1024,1,1024	1	1	+10:41:27.812	+82:56:33.423	0	0	0	0
5809	1.530918E9	1.530918E9	clear	ECMMD Amplifier	1	1	1,1024,1,1024	1	1	+11:05:25.176	+82:56:42.548	0	0	0	0
5810	1.530918E9	1.530918E9	clear	ECMMD Amplifier	1	1	1,1024,1,1024	1	1	+11:05:27.716	+82:56:42.561	0	0	0	0
5811	1.530919E9	1.530919E9	clear	ECMMD Amplifier	1	1	1,1024,1,1024	1	1	+20:45:54.867	+88:27:46.997	0	0	0	0
5812	1.530919E9	1.530919E9	clear	ECMMD Amplifier	1	1	1,1024,1,1024	1	1	+20:45:54.867	+88:27:46.997	0	0	0	0
5813	1.530919E9	1.530919E9	clear	ECMMD Amplifier	1	1	1,1024,1,1024	1	1	+20:46:23.701	+88:27:47.28	0	0	0	0
5814	1.530919E9	1.530919E9	clear	ECMMD Amplifier	1	1	1,1024,1,1024	1	1	+20:46:23.701	+88:27:47.28	0	0	0	0
5815	1.530920E9	1.530920E9	clear	ECMMD Amplifier	1	1	1,1024,1,1024	1	1	+20:54:28.998	+88:27:47.687	0	0	0	0
5816	1.530920E9	1.530920E9	clear	ECMMD Amplifier	1	1	1,1024,1,1024	1	1	+20:54:33.522	+88:27:47.694	0	0	0	0
5817	1.530920E9	1.530920E9	clear	ECMMD Amplifier	1	1	1,1024,1,1024	1	1	+21:00:38.237	+88:27:48.277	0	0	0	0
5818	1.530920E9	1.530920E9	clear	ECMMD Amplifier	1	1	1,1024,1,1024	1	1	+21:00:43.47	+88:27:48.295	0	0	0	0